

HIGH RESOLUTION CONVECTIVE HEAT TRANSFER MEASUREMENTS

Peter Ireland and Terry Jones

R-R UTC in Heat Transfer and
Aerodynamics

Dept. Engineering Science,
University of Oxford



Acknowledgements




- Internal cooling
 - David Gillespie, Calvin Tsang, Changmin Son
- Stepped solution and fin research
 - Andrew Neely
- Recuperator research underway
 - Marty Cerza and Juan Adams
- Transition work
 - Richard Anthony

Content

- High resolution htc measurements using temperature sensitive liquid crystals
 - Need for high resolution htc data
 - Scaling strategy
 - Liquid crystal instrument features
 - Application and test details
 - Example applications new developments
- Thin film gauges
 - Instrument details and recent developments
 - Applications
 - High density platinum gauges
- Conclusions

Need for high resolution htc data in turbomachinery



- Detailed thermal model of the engine component required for component life predictions
- Aerospace turbine blades are small and cooling systems are usually compact.



Example blade cooling temp distribution

- Blade far from isothermal
- Biot number not small enough

$$Biot = \frac{hL}{k}$$



Scaling issues




- Large scale model improves effective resolution.
- Switch off sideways (lateral) conduction to achieve local htc measurement with 1-d processing.
- No need for engine temperatures.
- Fluid dynamics correct through use of Reynolds number, Mach number and Prandtl number.

$$\text{Nu} = f(\text{Re}, \text{Pr}, \text{Mach})$$


Dimensionless heat transfer coefficient Dimensionless flow speed

The diagram illustrates the relationship between the Nusselt number (Nu) and its independent variables. The equation $\text{Nu} = f(\text{Re}, \text{Pr}, \text{Mach})$ is centered. Two arrows point from the labels "Dimensionless heat transfer coefficient" and "Dimensionless flow speed" towards the respective variables in the equation: "Nu" and "Mach".

Transient method with liquid crystals for internal cooling




Temperature measurement






Red Green Blue signals from video
→ Intensity or Hue processing

$$I = R + G + B$$

$$\cos(H) = \frac{2R - G - B}{\sqrt{6((R - I)^2 + (G - I)^2 + (B - I)^2)}}$$





Time of crystal colour change depends on local h



Heat transfer inside a film-cooling hole fed in cross-flow

Hole diameter = 0.3mm




Hole diameter = 22mm

:00:11:69


:00:16:19

:00:22:19

Recorded Colour Play within the Hole




Typical Intensity Histories at 6 Positions on the Film Cooling Hole Surface. Monochromatic processing.




htc

Local Nusselt Number Distribution, 90° Hole






Lattice cooling system for trailing edge

Geometric scaling essential
for resolution




The Assembled Lattice Model






Photographs of the Working Section




ONR Thermal Materials Workshop 2001


Nusselt number distribution at several Reynolds numbers






Re=17950 Nu(Dittus-Boelter) = 50.5




Re=25890 Nu(Dittus-Boelter)=67.7




Re=31060 Nu(Dittus-Boelter)=78.3





Re=17950 Nu(Dittus-Boelter)=50.5



Re=25890 Nu(Dittus-Boelter)=67.7




Re=31060 Nu(Dittus-Boelter)=78.3




Rib roughened passage

- 3 metres (60d) long, square cross section cooling passage
- Perspex walls with temperature sensitive liquid crystal coated on the inner surface.
- Reynolds number from 20,000 to 60,000
- Air is heated at the inlet using heater mesh
- Test section situated from 40d to 60d
- Fully developed flow established ahead of the test section
- Rotation not simulated in the experiment




Heat transfer distribution

60° interrupted inline ribs




Impingement heat transfer rig



The perspex test rig is instrumented with liquid crystal coated impingement and target plates, and fast response gas thermocouples at the entrance and exit of the working section

Impinging jet heat transfer




$$Re_{javg} = 26710$$




Data through the mesh heater

Reynolds Number = 38000

Image at time = 60 seconds



60 μm square apertures




Heater convective efficiency

- High convective efficiency (~50%) so:

- wires run cool
- radiation insignificant
- support easy to engineer

$$\eta_{CONVECTIVE} = \frac{T_{wire} - T_{downstream}}{T_{wire} - T_{upstream}}$$

- Suitable for switching temperature of low speed flows



Fluid temperature

Step change to T_g

Surface temperature

$$T_s = T_0 + (T_g - T_0) \left(1 - \exp\left(\frac{h^2 t}{\rho c k}\right) \times \operatorname{erfc}\left(\frac{h\sqrt{t}}{\sqrt{\rho c k}}\right) \right)$$

Series of steps

$$T_s = T_0 + \sum_{i=1}^n (T_{g_i} - T_{g_{i-1}}) \left(1 - \exp\left(\frac{h^2 (t - \tau_i)}{\rho c k}\right) \operatorname{erfc}\left(\frac{h\sqrt{t - \tau_i}}{\sqrt{\rho c k}}\right) \right)$$


Exponential with asymptote T_g

$$\begin{aligned} \frac{T_s - T_0}{T_g - T_0} &= I - \frac{\frac{\rho c k}{h^2 \tau}}{\left(I + \frac{\rho c k}{h^2 \tau} \right)} e^{-\frac{h^2 t}{\rho c k}} \\ \operatorname{erfc}\left(\frac{h\sqrt{t}}{\sqrt{\rho c k}}\right) &- e^{-\frac{t}{\tau}} \frac{I}{\left(I + \frac{\rho c k}{h^2 \tau} \right)} \\ &\left(I + \frac{\sqrt{\rho c k}}{h\sqrt{\tau}} \left(\frac{1}{\pi} \sqrt{\frac{t}{\tau}} + \frac{2}{\pi} \sum_{n=1}^{\infty} \frac{1}{n} e^{-\frac{n^2}{4}} \right) \right) \\ &\quad \sinh\left(n\sqrt{\frac{t}{\tau}}\right) \end{aligned}$$





Ramp function with slope m

$$T_s = T_0 + mt \left\{ 1 - \frac{2}{\beta} + \frac{1 - \exp(\beta^2) \operatorname{erfd}(\beta)}{\beta^2} \right\}$$







Works heating or cooling




Use of stepped flow temperature




Increasing the number of steps



Work on Dry Low Emission Combustor




Film Cooling Effectiveness



Lateral Conduction Correction

Step Two: Alternating Direction Methods (after Jim Douglas)

$$\frac{\partial T}{\partial t} = \alpha \left(\frac{\partial^2 T}{\partial x^2} + \frac{\partial^2 T}{\partial y^2} + \frac{\partial^2 T}{\partial z^2} \right)$$





- a) $(\Delta_x^2 - \frac{2}{\Delta t})T_{n+1}^* = -(\Delta_x^2 + 2\Delta_y^2 + 2\Delta_z^2 + \frac{2}{\Delta t})T_n$
- b) $(\Delta_y^2 - \frac{2}{\Delta t})T_{n+1}^{**} = \Delta_y^2 T_n - \frac{2}{\Delta t} T_{n+1}^*$
- c) $(\Delta_z^2 - \frac{2}{\Delta t})T_{n+1} = \Delta_z^2 T_n - \frac{2}{\Delta t} T_{n+1}^{**}$

where $\Delta_x^2 T_{i,j,k,n} = (T_{i+1,j,k,n} - 2T_{i,j,k,n} + T_{i-1,j,k,n}) / (\Delta x)^2$


Lateral Conduction

Heat Flux at each pixel can be then be calculated. The Film Cooling Effectiveness and Heat Transfer Coefficient can be determined by the interception of the

Temperature axis and the gradient of the heat flux graph shown above.



Current marine application



- 832 channels in the specimen.
- Length of each channel: 197 mm.

RR spiral recuperator




FIGURE 2




FIGURE 5






FIGURE 6



Earlier interrupted fin htc



Recuperator heat transfer research underway



Early thin film gauges



DIFFUSION
EQUATION


$$\frac{\partial^2 T}{\partial x^2} = \frac{\rho c}{k} \frac{\partial T}{\partial t}$$

LAPLACE
TRANSFORM

$$\dot{q} = \sqrt{\rho c k} \sqrt{sT}$$

FREQUENCY

$$\dot{q}(\omega) = \sqrt{\rho c k} \sqrt{j\omega T(\omega)}$$




Analysis

$$T \rightarrow \dot{q}$$

$$\bar{\dot{q}} = \sqrt{\rho c k} \sqrt{s T}$$

STEP $\dot{q}(t)$




$$\bar{\dot{q}} \sim \frac{1}{s}$$

$$\bar{T} \sim \frac{1}{s^{3/2}}$$

$$T(t) \sim t^{1/2}$$

ANALOGUE




$$\frac{\partial^2 v}{\partial x^2} = \frac{r}{c} \frac{\partial v}{\partial t}$$

$$V \equiv T$$

$$I \equiv \dot{q}$$

Two layer gauges





$$\bar{q} = \sqrt{\rho_1 c_1 k_1} \sqrt{sT}$$

$$x \left[\frac{1 - A \exp\left(-2l_1\left(\frac{s}{\alpha_1}\right)^{1/2}\right)}{1 + A \exp\left(-2l_1\left(\frac{s}{\alpha_1}\right)^{1/2}\right)} \right]$$





$$A = \frac{1-\sigma}{1+\sigma}$$


$$\sigma^2 = \frac{\rho_2 c_2 k_2}{\rho_1 c_1 k_1}$$




Gauge sensor (Platinum)
Copper leads C


Upilex
Glue






ONR Thermal Materials Workshop 2001







High Density Layout using Gauge Arrays




Model covered with sheet of TFG arrays



Turbulent Spot Detection with Surface TFG's



Visualizing Transitional Heat Flux



Conclusions

- Presentation of existing techniques for miniature measurements.
- Scaling can be useful for determination of convective loads at high resolution.
- Miniature thin film gauges being continuously developed.